

Human Proteins, Human Source[™]

Post-translational modifications on human cell expressed (HCXTM) human recombinant proteins

Wilson NL, Jiang H, Meier A, <u>Liddell K</u>

Apollo Cytokine Research Email: info@apollocytokineresearch.com

hox[™] Human Cell Expressed

www.apollocytokineresearch.com

INTRODUCTION

Most proteins undergo post-translational modification (PTM), which can alter their physical and chemical properties (e.g., MW, pl, folding, stability, activity, antigenicity, and function). The presence or absence of PTMs may be significant to both the activity and longevity of the protein in a biological system.

Various methods were used to study the differences in PTM, in particular glycosylation, between recombinant human proteins expressed in modified human 293 cells as opposed to non-human cells. These methods determine not only the differences in glycosylation but may also give some insight into the possible differences in function of the protein.

RESULTS AND DISCUSSION

1D-PAGE


Using enzymes which remove the oligosaccharides while leaving the protein backbone intact provides valuable information into the amount of glycosylation that occurs on a particular protein. This can be expressed as a percentage of the total mass of the protein as it appears on a onedimensional gel.


In the example shown here the glycoprotein appears as a fuzzy band with a molecular weight between 75-100 kDa. After removal of the N-linked oligosaccharides the protein resolves into a tighter band at 55 kDa, and subsequently drops to 53 kDa with the additional removal of the Olinked oligosaccharides. This indicates that approximately 25-45% of the glycoprotein mass is due to the N-linked oligosaccharides and 2% is due to the Olinked oligosaccharides.

CR013

Recombinant proteins are dependent on the machinery of the cell line in which they are made to determine their PTMs. Hence the PTMs of human proteins made recombinantly in a human cell line may differ significantly from the same protein made in NS0, CHO, E. coli or any other nonhuman cell line. For example *E. coli* does not possess the type of cellular machinery used for glycosylation in higher organisms, hence human proteins produced in an E. coli cell line are non-glycosylated. Consequently, the function of this protein may vary significantly from the glycosylated version.

One-dimensional electrophoresis combined with enzymatic deglycosylation was used to determine the relative mass of the glycosylated versus the deglycosylated protein. Enzymatic and chemical deglycosylation methods combined with MALDI-MS and LC-MS were used to determine the glycosylation sites as well as the N- and O-linked oligosaccharide structures present on the protein.





Lane 1: Prescision Plus Stds Lane 2: Glycosylated Protein Lane 3: De-N-glycosylated protein Lane 4: De-N&O-glycosylated protein Lane 5: Precision Plus Stds


1D and 2D-PAGE analysis of SCF sR

PEPTIDE MASS FINGERPRINTING

Human protein expressed in *E. coli* cell line


Human protein expressed in modified human 293 cell line

MATERIALS AND METHODS



(PMF) ANALYSIS

MS methods such as MALDI-MS typically are not able to detect glycopeptides due to their very large mass. By removing the oligosaccharides not only does the chance of seeing the peptides increase, but also critical information can be obtained as to the sites of glycosylation. Removing the Nlinked oligosaccharides using enzymatic methods changes the asparagine amino acid to an aspartic acid. This changes the mass of the peptide, hence enabling the determination of whether the site is non-, partially or completely glycosylated. Extra peptides visible after the removal of Olinked oligosaccharides, also enables the prediction of peptides where O-linked glycosylation may occur.






Table 1. PMF analysis from 1D gel

| Treatment | Tryptic Coverage | N-linked Sites | O-linked Sites |
|--|------------------|--------------------------------|---|
| none | 31% | 1 site seen not N-glycosylated | |
| PNGase F | 47% | 4 de-N-glycosylated sites seen | |
| PNGase F, Sialadase ß(1-4)Galactosidase | 39% | 4 de-N-glycosylated sites seen | 2 peptides not previously seen which contain 14 |

CONCLUSION


- Human cell expressed proteins differ other recombinant proteins from produced in non-human cell lines, such as NSO, E. coli and CHO, in the types of post-translational modifications that are present on the protein backbone.
- PTMs can alter the physical and chemical properties of the protein, and this may have a significant effect on the protein's activity, longevity, and immunogenicity in a biological system.
- E. coli does not possess the type of cellular machinery used for glycosylation, and hence human proteins produced in an E. coli cell line will be nonglycosylated.
- glycoproteomic Proteomic and 0 techniques can be used to determine the relative quantity, structure, and position of glycosylation, which help enhance our knowledge of how a particular protein may react in a biological system.

ß-N-Acetylglucosaminidase O-Glycanase

possible O-linked sites.

MS OLIGOSACCHARIDE STRUCTURAL ANALYSIS

Using the above methods to remove the oligosaccharides also enables the structure of the oligosaccharides to be determined using LC-MS/MS. This allows determination



of the oligosaccharide epitopes, which contribute to the function of the protein. Our studies have shown major differences in the glycosylation between human protein produced recombinantly from modified human 293 cells as opposed to CHO or NS0 cell lines.

